

COVID-19 Packaging & Shipping Using Kit 51

Prevent Disease – Promote Wellness – Improve Quality of Life



1

COVID-19 Kit 51 Components


- After kit receipt
 - Take ice pack out and place in freezer
 - Viral transport media can be stored at room temp. or refrigerated
- Note:** It's acceptable to keep kit at room temp except ice pack (freezer) until used



2

NP/OP Collection


- Use only synthetic fiber swabs with plastic shafts.
 - Do not use calcium alginate swabs or swabs with wooden shafts, as they may contain substances that inactivate some viruses and inhibit PCR testing.
- **Nasopharyngeal swab:** Insert a swab into the nostril parallel to the palate. Leave the swab in place for a few seconds to absorb secretions
- **Oropharyngeal swab** (e.g., throat swab): Swab the posterior pharynx, avoiding the tongue
- After collection, place swabs immediately into sterile tubes containing 2-3 ml of viral transport media (i.e. VTM, M₄).
 - NP and OP specimens may be kept in separate vials or **combined into a single vial.**



3

Packaging VTM Specimens

- Acceptable Specimens
 - NP Swabs & OP Swabs
- Packaging
 - Label each tube with: Patient name, Date of Birth, & Source i.e. NP/OP
 - **NOTE:** NP & OP swabs can **both** be placed in same VTM tube
 - Place VTM tube (tighten cap) in plastic bag with absorbent material square
 - Place sample(s) in 95kPa bag



4

Sputum Collection

- Have the patient rinse the mouth with water
- Next have patient expectorate deep cough sputum directly into a sterile, leak-proof, screw-cap sputum collection cup or sterile dry container.
 - Note: Sub-optimal samples will be rejected (if not from deep cough collection)
- Sputum (sputum can be collected in a variety of sterile containers
 - Sterile urine cup or 50ml conical tube
 - Note: Collection containers not provided in kit.

Sterile cup & conical tube



5

Packaging Sputum Specimens

- Acceptable Specimens
 - Sputum
- Packaging
 - Label container with:
 - Patient name, Date of Birth, & Source i.e. Sputum
 - Place sputum container inside the 95kPa bag with absorbent material square & VTM tubes.

Sputum, NP, & OP in separate bags with absorbent inside 95kPa bag



6

Complete all areas highlighted on the top portion of the form

DATE RECEIVED IN LABORATORY

LABORATORY SAMPLE NUMBER

Michigan Department of Health and Human Services - Bureau of Laboratories
P.O. Box 30035 3350 North Martin Luther King Jr. Blvd. Lansing, MI 48909
Laboratory Records: 517-335-8059 Technical Information: 517-335-8067
Fax: 517-335-9871 Web: www.michigan.gov/mdhhs/lab

SUBMITTER INFORMATION

AGENCY CODE (If known)

FP

TELEPHONE

FAX

NATIONAL PROVIDER IDENTIFIER #

CONTACT PERSON/ORDERING PHYSICIAN/PROVIDER NAME

PATIENT INFORMATION

NAME (LAST, FIRST, MIDDLE INITIAL) or UNIQUE IDENTIFIER

SUBMITTER PATIENT # (If Applicable)

CITY

ZIP

GENDER

RACE

ADAP NUMBER

BIRTH DATE (MM-DD-YYYY)

ETHNICITY

SUBMITTER SPECIMEN #

COLLECTION DATE (MM-DD-YY)

COLLECTION TIME (H:MM)

9

Indicate source in highlighted area

Check "Other" Box and write in: COVID-19

INDICATE TEST REQUESTED

INSTRUCTIONS FOR COMPLETION: Complete reverse side of form for corresponding numbers in parentheses and in bold.

INDICATE SPECIMEN SOURCE	SEROLOGY	MICROBIOLOGY	TESTS THAT REQUIRE MDHHS APPROVAL
☐ AMNIOTIC FLUID ☐ BRONCHIAL ☐ CERVIX ☐ CSF ☐ GASTRIC ☐ NASOPHARYNGEAL ☐ ORAL MUCOSAL TRANSUDATE ☐ PLASMA ☐ SERUM ☐ STOOL ☐ SPUTUM ☐ THROAT ☐ URETHRA ☐ URINE ☐ WHOLE BLOOD ☐ FOOD-Specific ☒ OTHER-Specific NP/OP Swab	SERUM STATUS - If Applicable ☐ ACUTE ☐ CONVALESCENT ☐ ARBOVIRUS ENCEP PANEL (IgM) May-Oct Includes: Eastern Equine, California, St. Louis and West Nile CSF Only ☐ BRUCELLE SEROLOGY ☐ FUNGAL SEROLOGY COMPLEMENT FIX ☐ FUNGAL IMMUNODIFFUSION ☐ FRANCISELLA SEROLOGY ☐ LEGIONELLA - HA ☐ LYME DISEASE - EIA (4) ☐ MEASLES IgG ☐ MUMPS IgG ☐ RABIES AB SEROLOGY (D) ☐ RUBELLA IgG ☐ TETANUS TOXIN EIA ☐ TETANUS TOXIN IgG ☐ SYPHILIS PANEL (1) ☐ SYPHILIS TP-PA (ONLY) (1) ☐ SYPHILIS VDRL - CSF Only (1) ☐ SYPHILIS DFA (1,2) ☐ SYPHILIS IgM WESTERN BLOT* (1)	☐ AEROBIC ISOLATE ID (5) ☐ ANTIMICROBIAL RESISTANCE CONF. (5) ☐ AFB SLIDE/CULTURE-CLINICAL SPECIMEN ☐ AFB IDENTIFICATION-ISOLATE ID ☐ ENTERIC BACTERIAL CULTURE ☐ FOODBORNE ILLNESS-Steal or Food (6) ☐ FUNGAL IDENTIFICATION - ISOLATE ID ☐ LEGIONELLA CULTURE ☐ NEISSERIA GONORRHOEA - ISOLATION ☐ NEISSERIA - REFERRED CULTURE ☐ PARASITOLOGY - BLOOD ☐ PARASITOLOGY - STOOL ☐ PARASITOLOGY - WOC ☐ PERTUSSIS PCR ☐ SALMONELLA SEROTYPING - HUMAN ☐ SHIGELLA SEROTYPING ☐ E. COLI SHIGA-TOXIN PRODUCER (STEC) ☐ ENTEROVIRUS PCR (6) ☐ RESPIRATORY PCR PANEL ☐ INFLUENZA (PCR/CULTURE) (7) PATIENT STATUS (Influenza) ☐ OUTPATIENT ☐ INPATIENT ☐ UNKNOWN ☐ VIRAL CULTURE	EMERGING ARBOVIRUS PANEL ☐ PCR ☐ SEROLOGY ☐ AFB NUCLEIC ACID AMPLIFICATION ☐ BACTERIAL TYPING-PFGE (6) ☐ BOTULISM TOXIN ☐ MUMPS - PCR ☐ MEASLES IgM ☐ NOROVIRUS PCR (6) ☐ PERTUSSIS CULTURE ☐ RUBELLA IgM (1) ☐ SALMONELLA SEROTYPING NON-HUMAN ☐ TOXIC SHOCK TESTING ☒ OTHER COVID-19 OTHER ☐ AUTOCLAVE TEST STRIPS ☐ LEGIONELLA - DFA ☐ LYME DISEASE - DFA (Tick) HEPATITIS TESTING ☐ HEPATITIS C ANTIBODY (1) ☐ HEPATITIS B SURFACE ANTIGEN (HbSAg) (1) ☐ HEPATITIS B ANTIBODY (Anti-HbSAg) (1) ☐ HEPATITIS A ANTIBODY (IgM) (1)

DCH-0853 (Rev. 01-19)

10

Packaging samples cont.

Place frozen ice pack
inside styrofoam
insert within the box



Place 95kPa bag with
samples on top of
frozen ice pack in box



11

Packaging samples cont.

- Place styrofoam lid on top
- Place completed State of Michigan-Laboratory Test Requisition on top of styrofoam lid

INDICATE TEST REQUESTED			
INSTRUCTIONS FOR COMPLETION: Complete reverse side of form for corresponding numbers in parentheses and in bold.			
INDICATE SPECIMEN SOURCE	SEROLOGY	MICROBIOLOGY	TESTS THAT REQUIRE MDHHS APPROVAL
☐ ANASTOTIC FLUID ☐ BRONCHIAL ☐ CERVIK ☐ CSF ☐ GASTRIC ☐ NASOPHARYNGEAL ☐ ORAL MUCOSAL TRANSDUCATE ☐ PLASMA ☐ SERUM ☐ STOOL ☐ SPUTUM ☐ THROAT ☐ URINE ☐ URINE ☐ WHOLE BLOOD ☒ NP/OP swab	☐ SERUM STATUS - F Agglutina ☐ ACUTE ☐ CONVALESCENT ☐ ARBOVIRUS ENCEP PANEL (IgM) ☐ May-Or Includes Eastern Equine Encephalitis ☐ IgG and IgM Encephalitis - CSF Only ☐ BRUGELLA SEROLOGY ☐ FUNGAL SEROLOGY COMPLEMENT FIX ☐ FUNGAL MANNOSE POLYMERIZATION ☐ FRANKIELLA SEROLOGY ☐ LEGIONELLA - IHA ☐ LYME DISEASE - EIA (4) ☐ MEASLES IgG ☐ MUMPS IgG ☐ RABIES AB SEROLOGY (S) ☐ RABIES IgG ☐ TETANUS TOXIN BA ☐ WARSZELLA ZOSTER IgG	☐ AEROBIC ISOLATE ID (S) ☐ ANTIMICROBIAL RESISTANCE COMP (S) ☐ AFB SLIDE CULTURE CLINICAL SPECIMEN ☐ AFB IDENTIFICATION ISOLATE ID ☐ ENTERIC BACTERIAL CULTURE ☐ FOODBORNE ILLNESS Swab w/ Food (S) ☐ FUNGAL IDENTIFICATION - ISOLATE ID ☐ LEGIONELLA CULTURE ☐ NEISSERIA GONORRHOEA - ISOLATION ☐ NEISSERIA - HEPERESIS CULTURE ☐ PARASITOLOGY - BLOOD ☐ PARASITOLOGY - STOOL ☐ PARASITOLOGY - WORM ☐ PERTUSSIS PCR ☐ PERTUSSIS PCR ☐ SALMONELLA SENSITIZING NON-HUMAN ☐ S. COOL (S) SENSITIZING PRODUCER (STEC)	☐ EMERGING ARBOVIRUS PANEL ☐ PCR ☐ SEROLOGY ☐ AFB NUCLEIC ACID AMPLIFICATION ☐ BACTERIAL TYPING PFGE (S) ☐ BOTULISM TOXIN ☐ HAVES - PCR ☐ MEASLES IgM ☐ MUMPS PCR (S) ☐ PERTUSSIS CULTURE ☐ RABIES IgM (S) ☐ SALMONELLA SENSITIZING NON-HUMAN ☒ EMERGING COVID-19 ☐ COVID-19 PCR TESTING OTHER ☐ AUTOLYSATE TEST STRIPS ☐ LEGIONELLA - OFA ☐ LYME DISEASE - OFA (T40) HEPATITIS TESTING ☐ HEPATITIS C ANTIBODY (S) ☐ HEPATITIS B SURFACE ANTIGEN ☐ HBsAg (S) ☐ HEPATITIS B ANTIBODY (Anti-HBc) (S) ☐ HEPATITIS B ANTIBODY (Anti-HBc) (S)
HIV TESTING ☐ HIV Ag/Ab - Serum (S) ☐ HIV Ag/Ab Oral Mucosal Transudate (S) ☐ CD4/CD8 (CD4 whole blood) (S) ☐ HIV-1 VIRAL LOAD (EDTA plasma) (S) ☐ HIV-1 GENOTYPING (EDTA plasma) (S)	SYPHILIS TESTING ☐ SYPHILIS PANEL (S) ☐ SYPHILIS TP-PA (QUALY) (S) ☐ SYPHILIS VDRL - CSF Only (S) ☐ SYPHILIS RPR (S) ☐ SYPHILIS IgG WESTERN BLOT (S)	VIROLOGY ☐ ENTEROVIRUS PCR (S) ☐ RESPIRATORY PCR PANEL ☐ INFLUENZA (PNCULTURE) (S) ☐ PATIENT STATUS (In-house) ☐ OUTPATIENT ☐ INPATIENT ☐ URINE ☐ URINE CULTURE	

12

Shipping Options

- Courier transport to BOL
 - Use Address label in kit # 51 with address:
3350 North Martin Luther King Jr. Blvd., Lansing
 - UPS (if courier is unavailable) Monday-Thursday only
 - Use UPS label included in kit # 51 with address:
927 Terminal Rd (our warehouse address)
 - UPS (if courier is unavailable) Friday only
 - Call the lab to have Friday overnight UPS label faxed to you
 - Contact (517) 335-8059
 - Weekend Delivery: Use your facility courier or refrigerate sample for Monday delivery
- Urgent requests call: Dr. Riner at (517) 230-7828



13

Shipping with Courier-Outer Box

- Close box and tape with packing tape.
- Place address label on top of sealed box
 - Fill in your name and facility address on label
- Place UN3373 Category B label on side of outer box
- Write full name and phone number with area code of Responsible Person from your facility on Top of box
- Place "refrigerate" yellow label on box



14

UPS Shipping-Monday-Friday ONLY

- Close box and tape with packing tape.
- Place UPS label on top of sealed box
- Place UN3373 Category B label on side of outer box
- Write full name and phone number with area code of Responsible Person from your facility on top of box
- Place “refrigerate” yellow label on box



15

Questions?

For questions regarding COVID-19 testing, contact:

Dr. Diana Riner

Virology Section Manager, MDHHS Bureau of Laboratories

rinerd@michigan.gov 517-335-8099

For questions for COVID-19 packaging:

Shannon Sharp

Bioterrorism Training Coordinator, MDHHS Bureau of Laboratories

sharpsh1@michigan.gov Office: 517-335-9653; cell (517) 331-7356

Matt Bashore

Supervisor DASH Unit, MDHHS Bureau of Laboratories

bashorem@michigan.gov Office: 517-335-8059;
Cell: 517-648-9804



16